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Now that you can measure a variety of properties of the SQUID in the Mr. SQUID probe, you can determine a key parameter of the device, namely the  $\beta_L$  parameter. This is defined by

$$\text{Eqn. 3-1} \quad \beta_L = \frac{2I_c L}{\Phi_0}$$

where  $L$  is the inductance of the SQUID loop.

Earlier we mentioned that the maximum modulation voltage depth  $\Delta V$  as measured from the  $V-\phi$  curve is related to the  $I_c R_N$  product that can be determined from the  $V-I$  curve. As we will discuss in Section 5, this relationship can be expressed simply in terms of the  $\beta_L$  parameter.

$$\text{Eqn. 3-2} \quad \beta_L = \frac{4I_c R_N}{\pi \Delta V} - 1$$

The above expression provides a simple way to determine the  $\beta_L$  parameter empirically (*i.e.*, without knowing the inductance  $L$  of the SQUID). This expression is strictly correct only if the critical currents of the two junctions are equal and only if thermal noise effects are negligible.

Both of these are approximations for the junctions in your Mr. SQUID.

The inductance of the Mr. SQUID chip may be written as the sum of four terms,

$L = L_h + L_k + L_{sl} + L_j$ , where  $L_h$  is the inductance of the SQUID washer,  $L_{sl}$  is the inductance of the long slit in the washer,  $L_k$  is the small kinetic inductance of the washer (arising from the

“inertia” of the electrons), and  $L_j$  is the inductance of the Josephson junction bridges, which also

includes a small kinetic inductance contribution. For a square washer with outer side length  $D$  and a square hole with inner side length  $d$  such that  $d/D < 1$ ,  $L_h = 1.24 \mu_0 d$ . Using  $d = 24 \mu\text{m}$  for the Mr. SQUID washer,  $L_h = 37.7 \text{ pH}$ . For the washer, the slit inductance per unit length is  $0.38 \text{ pH}/\mu\text{m}$  (it would be slightly less if the modulation coils were superconducting rather than normal metal). The slit in the Mr. SQUID washer is  $33 \mu\text{m}$  from the inner hole to the outer edge of the washer. Then,  $L_{sl} = 12.5 \text{ pH}$ . The kinetic inductance of the washer is more difficult to determine precisely but is estimated to be about  $2 \text{ pH}$ . The inductance of the Josephson junction bridges is estimated to be about  $8 \text{ pH}$ , which includes the kinetic inductance contribution (the inductance per unit length of the bridges is much higher than  $0.38 \text{ pH}/\mu\text{m}$  because of the narrow width of the bridges). Thus, the total inductance is approximately  $L = 60 \text{ pH}$ .

In general, for an  $N$ -turn coil integrated on top of the Mr. SQUID washer, the mutual inductance is given by

$$\text{Eqn. 3-3} \quad M = N [ L_h + L_k + 0.5 ( L_{sl} + L_j ) ]$$

For Mr. SQUID, the single-turn modulation coil actually goes around only three quarters of the washer (see Figure 3-6). Then, using  $N = 0.75$  as a rough approximation along with the inductance values given above,  $M = 37 \text{ pH}$ .

The horizontal axis of the  $V-\Phi$  curve measures the current through the modulation coil, and this is linearly related to the magnetic flux in the SQUID. The current gain in the Mr. SQUID control box

is 10,000 V/A (*i.e.*, 1 Volt =  $10^{-4}$  Amperes). The period of the modulation of the magnetic flux in the Mr. SQUID loop is the flux quantum,  $\Phi_0$  ( $2.07 \cdot 10^{-15}$  Wb in MKS units.)

By measuring the amount of current ( $\Delta I$ ) in the coil that is required to produce a change of one fluxon through the hole of the SQUID, you can determine the mutual inductance ( $M$ ) of the SQUID, using the following formula:

Eqn. 3-4 
$$M = \frac{\Phi_0}{\Delta I}$$

You can find the value of  $\beta_L$  for the SQUID in your Mr. SQUID probe using Eqn. you're your measurement of  $I_C$ , and  $L \approx 60$  pH. From the measured values of  $I_C$  and  $\Delta V$ , you can calculate  $\beta_L$  using Eqn. 3-2. Compare these two values. Do they agree?

The fact that the values calculated using Eqn. 3-1 and Eqn. 3-2 do not agree was a mystery for a number of years after the 1986 discovery of high- $T_C$  materials, as these equations worked quite well for SQUIDs made using traditional low- $T_C$  materials. This lack of agreement was resolved in 1993 with the recognition that thermal effects were playing a large role in the behavior of high- $T_C$  SQUIDs. The lack of agreement between Eqn. 3-1 and Eqn. 3-2 is in large part due to the fact that the Mr. SQUID is at a relatively high temperature where thermal energies ( $k_B T$ ) are no longer small compared to the energy of a flux quantum ( $\Phi_0^2/L$ ). The relationship between the observed voltage modulation and  $\beta_L$  at a nonzero temperature  $T$  changes Eqn. 3-2 to :

Eqn. 3-5 
$$\beta_L = \frac{4I_C R_N}{\pi \Delta V} \left( 1 - 3.57 \frac{\sqrt{k_B T L}}{\Phi_0} \right) - 1$$

Figure 3-9 illustrates the differences between the  $\beta_L$  values calculated using Eqn. 3-1, Eqn. 3-2, and Eqn. 3-5.


Figure 3-9 Values of  $\beta = 2I_C L/\Phi_0$  versus measurements of  $\beta_L$  based upon  $I_C R_N/\Delta V$  for 44 Mr.SQUID probes.

As one can see from Figure 3-9, the method of calculating  $\beta_L$  that takes thermal effects into account agrees quite well with the inductive  $\beta_L$  measurements.

On the following page, we include a data sheet for your SQUID. You may wish to use it (or better still, a photocopy) to enter your measurements on Mr. SQUID. An interesting and important issue for SQUID applications is how the SQUID sensor varies with age. You can track this with Mr. SQUID by filling out a copy of the data sheet each time you use the system.

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